

Equip. ID	Maker	Model	process	Serial No.
SWOS02	SES 200mm	Sugai Special	Solvent Strip Station	O-1996

Bath configurations: SWOS02							
Position No	1	2	3	4	5	6	7
Bath	CC	ORG	ORG	IPA	QDR	F/R	IPA DRYER

Summary of inspection results

Items inspected	Conditions of electric cable	Bath and bath parts	Parts (valve,sensor)	Operation Manual	Conditions of the units	Sub-module included
Results	All cables are OK (PLC also OK)	All Bath and parts confirmed	All Parts are OK (No missing parts)	Not found, but it might be inside of tool	No missing parts. Cable is OK. Plumbings also OK	

Packing List

SWOS02

SWOS02	Sugai Special (200mm)	Results inspected
Pack No.	Item	
1	CC, Organic Module	OK
2	IPA Module	OK
3	LD Module	OK
4	Main Electric Box	OK
5	F/R, IPA Dryer Module	OK
6	ULD Module	OK
7	Main PC Module	OK
8	IPA Dryer	OK
9	Circulation unit	OK
10	Org 1 circulation unit	OK
11	Org 2 circulation unit	OK
12	Conveyer # 1	OK
13	Blue Box (MISC, Etc)	OK
14	Monitor PC Rack	OK
15	Hepa filter unit	OK
16	Electric box	OK
17	IPA Scrubber	OK
18	Blue Box (I/D Control Box)	OK
19	Blue Box (Org 2 Quartz bath)	OK
20	Electric box 2	OK
21	Blue Box (Org 1 Quartz bath)	OK
22	Blue Box (IPA Quartz bath, tube)	OK
23	Conveyer # 2	OK
24	Blue Box (Quartz tube)	OK
25	Blue Box (Chuck 4ea)	OK
26	Conveyer # 3	OK
27	Conveyer # 4	OK
28	Quartz Dry (IPA)	OK
29	Hepa filter unit	OK
30	Blue Box (Dual power box)	OK



Bath configuration of SWOS02



Main body of SWOS02



Electric Rack and bath inside

